

How To Use Your Resound Vivia, Nexia, Savi or Enzo hearing aids at the Fox Theatre in Atlanta


The Fox Theatre in Atlanta has an Auracast™ audio system that allows you to walk into the venue and connect instantly to stream the performance audio directly through your hearing aids for crystal-clear sound.

Requirements:



1. Resound Vivia, Nexia, Enzo or Savi hearing aids
2. Your hearing aids must be paired with your smartphone using the Resound Smart 3D app.

Instructions:


1. Open your Resound Smart 3D app
2. Tap on “Find a broadcast”




3. Select the broadcast “Fox Theatre” from the list



4. Sound should start streaming directly through your hearing aids when you see the green icon on the right.



5. You can **adjust the volume of the Fox Theatre audio stream** on the top volume control. You can **adjust the volume of the Surroundings** which changes the volume of your hearing aid microphones on top of your ears (ex. someone behind you is eating popcorn or talking and you want to turn them down).



6. Please note that the broadcast will disconnect at intermission when you leave the theatre. You will need to reconnect when you return to your seat.

