




TIP: How to add remote control function

When you swipe up on your home screen, the  icon will allow you to make adjustments to your hearing aids.



When you swipe up on your home screen...



...this "Control Center" will pop up.

Follow these simple steps to set up this feature.

STEP 1: Tap **Settings**



on your home screen.


STEP 2: Under Settings, tap **Control Center**




STEP 3: Under Control Center, tap **Customize Controls**




STEP 4: Scroll down to “**MFi Hearing Devices**” and tap the 




This will move “MFi Hearing Devices” to the top portion.



STEP 5: When you swipe up, tap the functions on your hearing aids.



icon to access the remote control




Keep in mind this battery level indicator is not 100% accurate of your battery level and simply informs you the devices are connected.

Change the volume by swiping the bar left (down) or right (up).

Tap here to change your program. **(You may not have this option if your hearing aids are programmed to be entirely automatic.)**

If you have additional programs, they will appear like this:



NOTE: If you iPhone model has a “home” button, you can also click the button **THREE** times to access the remote control function.