



# HEARING CENTER


Eye Consultants of Atlanta

How to pair your hearing aids to your iPhone at home


## Step 1: Go to Settings App




## Step 2: Go to Accessibility




## Step 3: Go to Hearing Devices




**Step 4:** Turn off your hearing aids by pressing buttons on BOTH aids SIMULTANEOUSLY for 5 seconds, then turn BOTH back on by pressing BOTH buttons another 5 seconds. If aids are non-rechargeable, open & close battery doors instead of pushing buttons.



**Step 5:** Your name should appear. Once you see (R+L), then tap on your name to connect aids to the iPhone.




**Step 6:** Tap on Pair, once for each hearing aid




Turn page over



**Step 7:** Push "< hearing devices", now that aids are connected



**Step 8:** You will see it says both aids are connected



**TIP:** Turn off play ringtones so that when the phone rings, it rings from your phone and NOT in your hearing aids. This allows you to be able to find your phone and is less likely to frighten you!

**Tip:** Also, leave Hearing aid compatibility OFF!

